




COLOR RENDITION TEST CHART

REFLECTANCE



The test chart consists of

- a neutral gray background of 18% remission
- a white field of 89.9% remission
- a black field (Super Black Hole) with a remission of $< 0.5\%$
- a second black field, where color patterns may be applied alternatively

For the color pattern we supply 15 small plastic plates which are covered with self-adhesive foil on one side. Individual color patterns can be put easily on these plates. The self-adhesive foil is designed for manual use and its structure allows blisters that may occur to be removed by smoothing out. We recommend to cover the color pattern with a protective paper. The color pattern should be mounted on only that part of the plate, where the foil is. Colormetric measurement should be effected after mounting. The plastic plates with color patterns can be fixed on the test chart with the clamp.